

Galectin-8

Catalog No: #AP85138



Package Size: #AP85138-1 50ug #AP85138-2 100ug #AP85138-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galectin-8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Gal-8,Po66 carbohydrate-binding protein,Po66-CBP,Prostate carcinoma tumor antigen 1,PCTA-1
Accession No.	Uniprot ID:O00214
Uniprot	O00214
GeneID	3964;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LGALS8

Note: This product is for in vitro research use only