

Leukemia inhibitory factor receptor

Catalog No: #AP85140



Package Size: #AP85140-1 50ug #AP85140-2 100ug #AP85140-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Leukemia inhibitory factor receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	LIF receptor,CD118
Accession No.	Uniprot ID:O70535
Uniprot	O70535
GeneID	81680;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lifr

Note: This product is for in vitro research use only