

Protein max

Catalog No: #AP85141



Package Size: #AP85141-1 50ug #AP85141-2 100ug #AP85141-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein max
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Myc-associated factor X, Myc-binding novel HLH/LZ protein, Protein myn, Myn
Accession No.	Uniprot ID:P28574
Uniprot	P28574
GeneID	17187;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Max

Note: This product is for in vitro research use only