

# Phosphoribosyl pyrophosphate synthase-associated protein 1

Catalog No: #AP85144

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP85144-1 50ug #AP85144-2 100ug #AP85144-3 1mg

## Description

Product Name	Phosphoribosyl pyrophosphate synthase-associated protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PRPP synthase-associated protein 1
Accession No.	Uniprot ID:Q08DW2
Uniprot	Q08DW2
GeneID	789409;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PRPSAP1

Note: This product is for in vitro research use only