

Galectin-9

Catalog No: #AP85145



Package Size: #AP85145-1 50ug #AP85145-2 100ug #AP85145-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Galectin-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Gal-9
Accession No.	Uniprot ID:O08573
Uniprot	O08573
GeneID	16859;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Lgals9

Note: This product is for in vitro research use only