

## Creatine kinase M-type

Catalog No: #AP85146



Package Size: #AP85146-1 50ug #AP85146-2 100ug #AP85146-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Creatine kinase M-type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Creatine kinase M chain,M-CK
Accession No.	Uniprot ID:Q5XLD3
Uniprot	Q5XLD3
GeneID	397264;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CKM

Note: This product is for in vitro research use only