

# Methionine synthase

Catalog No: #AP85148



Package Size: #AP85148-1 50ug #AP85148-2 100ug #AP85148-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Methionine synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	5-methyltetrahydrofolate--homocysteine methyltransferase,Vitamin-B12 dependent methionine synthase,MS
Accession No.	Uniprot ID:A6H5Y3
Uniprot	A6H5Y3
GeneID	238505;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Mtr

Note: This product is for in vitro research use only