

VASP (Phospho-Thr278) Antibody

Catalog No: #12699



Package Size: #12699-1 50ul #12699-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VASP (Phospho-Thr278) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	VASP (Phospho-Thr278) Antibody detects endogenous levels of VASP only when phosphorylated at Thr278
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human VASP (Phospho-Thr278)
Target Name	VASP
Modification	Phospho
Other Names	VASP
Accession No.	Swiss-Prot#: P50552NCBI Gene ID: 7408
Uniprot	P50552
GeneID	7408;
Target Species	human
Calculated MW	46kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only