

## Matrix metalloproteinase-14

Catalog No: #AP85150



Package Size: #AP85150-1 50ug #AP85150-2 100ug #AP85150-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Matrix metalloproteinase-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	MMP-14, Membrane-type matrix metalloproteinase 1, MT-MMP 1, Membrane-type-1 matrix metalloproteinase, MT1-MMP
Accession No.	Uniprot ID:Q9GLE4
Uniprot	Q9GLE4
GeneID	281915;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: MMP14

Note: This product is for in vitro research use only