

Multidrug resistance-associated protein 1

Catalog No: #AP85152



Package Size: #AP85152-1 50ug #AP85152-2 100ug #AP85152-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Multidrug resistance-associated protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ATP-binding cassette sub-family C member 1,Leukotriene C(4) transporter,LTC4 transporter,MRP1
Accession No.	Uniprot ID:Q8HXQ5
Uniprot	Q8HXQ5
GeneID	281588;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ABCC1

Note: This product is for in vitro research use only