

Beta-microseminoprotein

Catalog No: #AP85153



Package Size: #AP85153-1 50ug #AP85153-2 100ug #AP85153-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-microseminoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Prostate secreted seminal plasma protein, Prostate secretory protein of 94 amino acids, PSP-94, Psp94
Accession No.	Uniprot ID:P97580
Uniprot	P97580
GeneID	29311;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Msmb

Note: This product is for in vitro research use only