

Multiple coagulation factor deficiency protein 2 homolog

Catalog No: #AP85156

Package Size: #AP85156-1 50ug #AP85156-2 100ug #AP85156-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Multiple coagulation factor deficiency protein 2 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Neural stem cell-derived neuronal survival protein,Sdnst
Accession No.	Uniprot ID:Q8K5B2
Uniprot	Q8K5B2
GeneID	193813;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mcf2

Note: This product is for in vitro research use only