

Androgen Receptor (Phospho-Ser308) Antibody

Catalog No: #12700



Package Size: #12700-1 50ul #12700-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Androgen Receptor (Phospho-Ser308) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	Androgen Receptor (Phospho-Ser308) Antibody detects endogenous levels of Androgen Receptor only when phosphorylated at AnSer308
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Androgen Receptor (Phospho-Ser308)
Target Name	Androgen Receptor
Modification	Phospho
Other Names	AR, Androgen receptor, Dihydrotestosterone receptor, DHTR, AIS, KD, NR3C4, HUMARA, SBMA, TFM, SMAX1, HYSP1
Accession No.	Swiss-Prot#: P10275NCBI Gene ID: 367
Uniprot	P10275
GeneID	367;
Target Species	human
Calculated MW	98kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of Hela cell (Left) and Hela cell Serum-free treated (Right), using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only