

La-related protein 7

Catalog No: #AP85160

Package Size: #AP85160-1 50ug #AP85160-2 100ug #AP85160-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	La-related protein 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	La ribonucleoprotein domain family member 7,D3Wsu161e
Accession No.	Uniprot ID:Q05CL8
Uniprot	Q05CL8
GeneID	28036;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Larp7

Note: This product is for in vitro research use only