

Lysyl oxidase homolog 3

Catalog No: #AP85165



Package Size: #AP85165-1 50ug #AP85165-2 100ug #AP85165-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lysyl oxidase homolog 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Lysyl oxidase-like protein 3,LOXL
Accession No.	Uniprot ID:P58215
Uniprot	P58215
GeneID	84695;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LOXL3

Note: This product is for in vitro research use only