

Lutropin-choriogonadotropic hormone receptor

Catalog No: #AP85166



Package Size: #AP85166-1 50ug #AP85166-2 100ug #AP85166-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lutropin-choriogonadotropic hormone receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	LH/CG-R,Luteinizing hormone receptor,LSH-R
Accession No.	Uniprot ID:Q90674
Uniprot	Q90674
GeneID	395776;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LHCGR

Note: This product is for in vitro research use only