

Interferon-stimulated gene 20 kDa protein

Catalog No: #AP85170



Package Size: #AP85170-1 50ug #AP85170-2 100ug #AP85170-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon-stimulated gene 20 kDa protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Estrogen-regulated transcript 45 protein, Promyelocytic leukemia nuclear body-associated protein ISG20, HEM45
Accession No.	Uniprot ID: Q96AZ6
Uniprot	Q96AZ6
GeneID	3669;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ISG20

Note: This product is for in vitro research use only