

Protein Mdm4

Catalog No: #AP85171



Package Size: #AP85171-1 50ug #AP85171-2 100ug #AP85171-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein Mdm4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Double minute 4 protein,Mdm2-like p53-binding protein,Protein Mdmx,p53-binding protein Mdm4,MDMX
Accession No.	Uniprot ID:O15151
Uniprot	O15151
GeneID	4194;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MDM4

Note: This product is for in vitro research use only