

## Meprin A subunit alpha

Catalog No: #AP85172



Package Size: #AP85172-1 50ug #AP85172-2 100ug #AP85172-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Meprin A subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Endopeptidase-2,Endopeptidase-24.18 subunit alpha,E-24.18,MEP-1
Accession No.	Uniprot ID:Q64230
Uniprot	Q64230
GeneID	25684;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Mep1a

Note: This product is for in vitro research use only