

Lysyl oxidase homolog 2

Catalog No: #AP85174



Package Size: #AP85174-1 50ug #AP85174-2 100ug #AP85174-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysyl oxidase homolog 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Lysyl oxidase-like protein 2,Lysyl oxidase-related protein 2,Lysyl oxidase-related protein WS9-14
Accession No.	Uniprot ID:Q9Y4K0
Uniprot	Q9Y4K0
GeneID	4017;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LOXL2

Note: This product is for in vitro research use only