

Membrane cofactor protein

Catalog No: #AP85176



Package Size: #AP85176-1 50ug #AP85176-2 100ug #AP85176-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Membrane cofactor protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TLX, Trophoblast leukocyte common antigen, CD46, MCP, MIC10
Accession No.	Uniprot ID:P15529
Uniprot	P15529
GeneID	4179;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CD46

Note: This product is for in vitro research use only