

Lysosomal alpha-glucosidase

Catalog No: #AP85177



Package Size: #AP85177-1 50ug #AP85177-2 100ug #AP85177-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lysosomal alpha-glucosidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Acid maltase
Accession No.	Uniprot ID:Q9MYM4
Uniprot	Q9MYM4
GeneID	280798;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GAA

Note: This product is for in vitro research use only