

E-selectin

Catalog No: #AP85178



Package Size: #AP85178-1 50ug #AP85178-2 100ug #AP85178-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E-selectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CD62 antigen-like family member E,Endothelial leukocyte adhesion molecule 1,ELAM-1,Leukocyte-endothelial cell adhesion molecule 2,LECAM2,CD62E
Accession No.	Uniprot ID:P98110
Uniprot	P98110
GeneID	397508;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SELE

Note: This product is for in vitro research use only