

Mitogen-activated protein kinase 14

Catalog No: #AP85188



Package Size: #AP85188-1 50ug #AP85188-2 100ug #AP85188-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase 14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAP kinase 14, Cytokine suppressive anti-inflammatory drug-binding protein, CSAID-binding protein, MAP kinase MXI2, MAX-interacting protein 2, Mitogen-activated protein kinase p38 alpha, MAP kinase p38 alpha, Stress-activated protein kinase 2a, SAPK2a, CSBP, CSBP1, CSBP2, CSPB1, MXI2, SAPK2A
Accession No.	Uniprot ID:Q16539
Uniprot	Q16539
GeneID	1432;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MAPK14

Note: This product is for in vitro research use only