

Lysosomal-associated transmembrane protein 4B

Catalog No: #AP85191



Package Size: #AP85191-1 50ug #AP85191-2 100ug #AP85191-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Lysosomal-associated transmembrane protein 4B |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Accession No. | Uniprot ID:Q91XQ6 |
| Uniprot | Q91XQ6 |
| GeneID | 114128; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Laptm4b

Note: This product is for in vitro research use only