

## L-lactate dehydrogenase A chain

Catalog No: #AP85192



Package Size: #AP85192-1 50ug #AP85192-2 100ug #AP85192-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	L-lactate dehydrogenase A chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	LDH-A,LDH muscle subunit,LDH-M,Ldh-1,Ldh1
Accession No.	Uniprot ID:P04642
Uniprot	P04642
GeneID	24533;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ldha

Note: This product is for in vitro research use only