

## Serine protease inhibitor Kazal-type 7

Catalog No: #AP85196



Package Size: #AP85196-1 50ug #AP85196-2 100ug #AP85196-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine protease inhibitor Kazal-type 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Esophagus cancer-related gene 2 protein, ECRG-2, Ecg2
Accession No.	Uniprot ID:Q6IE32
Uniprot	Q6IE32
GeneID	408198;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Spink7

Note: This product is for in vitro research use only