

Melanocortin receptor 5

Catalog No: #AP85197



Package Size: #AP85197-1 50ug #AP85197-2 100ug #AP85197-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanocortin receptor 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MC5-R,MC-2
Accession No.	Uniprot ID:P33032
Uniprot	P33032
GeneID	4161;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MC5R

Note: This product is for in vitro research use only