

## Lithostathine-1

Catalog No: #AP85198



Package Size: #AP85198-1 50ug #AP85198-2 100ug #AP85198-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lithostathine-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Islet of Langerhans regenerating protein 1,REG 1,Pancreatic stone protein 1,PSP,Pancreatic thread protein 1,PTP,Regenerating protein 1
Accession No.	Uniprot ID:P43137
Uniprot	P43137
GeneID	19692;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Reg1

Note: This product is for in vitro research use only