

## Protein ERGIC-53

Catalog No: #AP85199



Package Size: #AP85199-1 50ug #AP85199-2 100ug #AP85199-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein ERGIC-53
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ER-Golgi intermediate compartment 53 kDa protein, Lectin mannose-binding 1, p58, Ergic53
Accession No.	Uniprot ID:Q9D0F3
Uniprot	Q9D0F3
GeneID	70361;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Lman1

Note: This product is for in vitro research use only