

Mitofusin-1

Catalog No: #AP85200



Package Size: #AP85200-1 50ug #AP85200-2 100ug #AP85200-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mitofusin-1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Mitochondrial transmembrane GTPase FZO1B, Transmembrane GTPase MFN1, Fzo1b |
| Accession No. | Uniprot ID:Q8R4Z9 |
| Uniprot | Q8R4Z9 |
| GeneID | 192647; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mfn1

Note: This product is for in vitro research use only