

Major facilitator superfamily domain-containing protein 2A



Catalog No: #AP85201

Orders: order@signalwayantibody.com

Package Size: #AP85201-1 50ug #AP85201-2 100ug #AP85201-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Major facilitator superfamily domain-containing protein 2A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mfsd2
Accession No.	Uniprot ID:Q9DA75
Uniprot	Q9DA75
GeneID	76574;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mfsd2a

Note: This product is for in vitro research use only