

Maternal embryonic leucine zipper kinase

Catalog No: #AP85202



Package Size: #AP85202-1 50ug #AP85202-2 100ug #AP85202-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Maternal embryonic leucine zipper kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein kinase PK38,mPK38,Tyrosine-protein kinase MELK,Kiaa0175,Pk38
Accession No.	Uniprot ID:Q61846
Uniprot	Q61846
GeneID	17279;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Melk

Note: This product is for in vitro research use only