

Phosphatidylcholine translocator ABCB4

Catalog No: #AP85203



Package Size: #AP85203-1 50ug #AP85203-2 100ug #AP85203-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphatidylcholine translocator ABCB4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ATP-binding cassette sub-family B member 4, Multidrug resistance protein 2, Multidrug resistance protein 3, P-glycoprotein 2, P-glycoprotein 3, Mdr2, Pgp3, Pgy2
Accession No.	Uniprot ID:Q08201
Uniprot	Q08201
GeneID	24891;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Abcb4

Note: This product is for in vitro research use only