

N-lysine methyltransferase KMT5A

Catalog No: #AP85204



Package Size: #AP85204-1 50ug #AP85204-2 100ug #AP85204-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-lysine methyltransferase KMT5A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	H4-K20-HMTase KMT5A, Histone-lysine N-methyltransferase KMT5A, Lysine N-methyltransferase 5A, Lysine-specific methylase 5A, PR/SET domain-containing protein 07, PR-Set7, SET domain-containing protein 8, PRSET7, SET07, SET8, SETD8
Accession No.	Uniprot ID:Q9NQR1
Uniprot	Q9NQR1
GeneID	387893;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: KMT5A

Note: This product is for in vitro research use only