

Laminin subunit alpha-5

Catalog No: #AP85205



Package Size: #AP85205-1 50ug #AP85205-2 100ug #AP85205-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Laminin subunit alpha-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Laminin-10 subunit alpha,Laminin-11 subunit alpha,Laminin-15 subunit alpha
Accession No.	Uniprot ID:Q61001
Uniprot	Q61001
GeneID	16776;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lama5

Note: This product is for in vitro research use only