

Neutrophil collagenase

Catalog No: #AP85206



Package Size: #AP85206-1 50ug #AP85206-2 100ug #AP85206-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Neutrophil collagenase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Matrix metalloproteinase-8,MMP-8 |
| Accession No. | Uniprot ID:O88766 |
| Uniprot | O88766 |
| GeneID | 63849; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mmp8

Note: This product is for in vitro research use only