

Dual specificity mitogen-activated protein kinase kinase 7

Catalog No: #AP85207

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP85207-1 50ug #AP85207-2 100ug #AP85207-3 1mg

Description

Product Name	Dual specificity mitogen-activated protein kinase kinase 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAP kinase kinase 7, JNK-activating kinase 2, MAPK/ERK kinase 7, MEK 7, Stress-activated protein kinase kinase 4, SAPK kinase 4, c-Jun N-terminal kinase kinase 2, JNK kinase 2, JNKK2, MEK7, MKK7, PRKMK7, SKK4
Accession No.	Uniprot ID:O14733
Uniprot	O14733
GeneID	5609;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MAP2K7

Note: This product is for in vitro research use only