

Metallo-beta-lactamase domain-containing protein 1

Catalog No: #AP85209



Package Size: #AP85209-1 50ug #AP85209-2 100ug #AP85209-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Metallo-beta-lactamase domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q2HJB0
Uniprot	Q2HJB0
GeneID	505636;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MBLAC1

Note: This product is for in vitro research use only