

Mitogen-activated protein kinase 15

Catalog No: #AP85213

Package Size: #AP85213-1 50ug #AP85213-2 100ug #AP85213-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase 15
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MAP kinase 15, Extracellular signal-regulated kinase 7, ERK-7, Erk7
Accession No.	Uniprot ID:Q9Z2A6
Uniprot	Q9Z2A6
GeneID	286997;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mapk15

Note: This product is for in vitro research use only