

E-selectin

Catalog No: #AP85215



Package Size: #AP85215-1 50ug #AP85215-2 100ug #AP85215-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E-selectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CD62 antigen-like family member E,Endothelial leukocyte adhesion molecule 1,ELAM-1,Leukocyte-endothelial cell adhesion molecule 2,LECAM2,CD62E,Elam-1
Accession No.	Uniprot ID:P98105
Uniprot	P98105
GeneID	25544;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sele

Note: This product is for in vitro research use only