

Leucine-rich repeat-containing G-protein coupled receptor 4

Catalog No: #AP85217

Package Size: #AP85217-1 50ug #AP85217-2 100ug #AP85217-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Leucine-rich repeat-containing G-protein coupled receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	G-protein coupled receptor 48,Gpr48
Accession No.	Uniprot ID:Q9Z2H4
Uniprot	Q9Z2H4
GeneID	286994;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lgr4

Note: This product is for in vitro research use only