

Vesicular integral-membrane protein VIP36

Catalog No: #AP85218



Package Size: #AP85218-1 50ug #AP85218-2 100ug #AP85218-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vesicular integral-membrane protein VIP36
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Lectin mannose-binding 2,Vesicular integral-membrane protein 36,VIP36
Accession No.	Uniprot ID:P49256
Uniprot	P49256
GeneID	403938;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LMAN2

Note: This product is for in vitro research use only