

Arachidonate 15-lipoxygenase

Catalog No: #AP85220



Package Size: #AP85220-1 50ug #AP85220-2 100ug #AP85220-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Arachidonate 15-lipoxygenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	15-LOX,12/15-lipoxygenase,Arachidonate 12-lipoxygenase, leukocyte-type,12-LOX,Arachidonate omega-6 lipoxygenase,LOG15
Accession No.	Uniprot ID:P16050
Uniprot	P16050
GeneID	246;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ALOX15

Note: This product is for in vitro research use only