

## Lipid phosphate phosphohydrolase 1

Catalog No: #AP85221



Package Size: #AP85221-1 50ug #AP85221-2 100ug #AP85221-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lipid phosphate phosphohydrolase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	35 kDa PAP,mPAP,Hydrogen peroxide-inducible protein 53,Hic53,PAP2-alpha,Phosphatidate phosphohydrolase type 2a,Phosphatidic acid phosphatase 2a,PAP-2a,PAP2a,Hpic53,Lpp1
Accession No.	Uniprot ID:Q61469
Uniprot	Q61469
GeneID	19012;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ppap2a

Note: This product is for in vitro research use only