

# Lymphocyte antigen 6 complex locus protein G6f

Catalog No: #AP85222



Package Size: #AP85222-1 50ug #AP85222-2 100ug #AP85222-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lymphocyte antigen 6 complex locus protein G6f
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	G6F
Accession No.	Uniprot ID:Q0V881
Uniprot	Q0V881
GeneID	539920;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LY6G6F

Note: This product is for in vitro research use only