

Membrane cofactor protein

Catalog No: #AP85224



Package Size: #AP85224-1 50ug #AP85224-2 100ug #AP85224-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Membrane cofactor protein |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Pig |
| Other Names | CD46,MCP |
| Accession No. | Uniprot ID:O02839 |
| Uniprot | O02839 |
| GeneID | 396922; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD46

Note: This product is for in vitro research use only