

## DAPK1 (Phospho-Ser308) Antibody

Catalog No: #12707



Package Size: #12707-1 50ul #12707-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DAPK1 (Phospho-Ser308) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms
Specificity	DAPK1 (Phospho-Ser308) Antibody detects endogenous levels of DAPK1 only when phosphorylated at Ser308
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human DAPK1 (Phospho-Ser308)
Target Name	DAPK1
Modification	Phospho
Other Names	DAPK1, DAP kinase 1, DAPK
Accession No.	Swiss-Prot#: P53355NCBI Gene ID: 1612
Uniprot	P53355
GeneID	1612;
Target Species	human
Calculated MW	160kd
Concentration	1.0mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

## Application Details

WB (1/500 - 1/1000)

## Images



Western blot analysis of DAP Kinase 1 (pS308) expression in H1792 (A), A549 (B) whole cell lysates.

---

Note: This product is for in vitro research use only