

# Methionine synthase

Catalog No: #AP85230



Package Size: #AP85230-1 50ug #AP85230-2 100ug #AP85230-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                                                         |
|--------------------|---------------------------------------------------------------------------------------------------------|
| Product Name       | Methionine synthase                                                                                     |
| Brief Description  | Recombinant Protein                                                                                     |
| Host Species       | E.coli                                                                                                  |
| Purification       | Greater than 90% by SDS-PAGE                                                                            |
| Species Reactivity | Rt                                                                                                      |
| Other Names        | 5-methyltetrahydrofolate--homocysteine methyltransferase, Vitamin-B12 dependent methionine synthase, MS |
| Accession No.      | Uniprot ID:Q9Z2Q4                                                                                       |
| Uniprot            | Q9Z2Q4                                                                                                  |
| GeneID             | 81522;                                                                                                  |
| Tag Info           | His                                                                                                     |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                                                  |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Mtr

Note: This product is for in vitro research use only