

Hepatic triacylglycerol lipase

Catalog No: #AP85232



Package Size: #AP85232-1 50ug #AP85232-2 100ug #AP85232-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hepatic triacylglycerol lipase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HL,Lipase member C,HTGL
Accession No.	Uniprot ID:P11150
Uniprot	P11150
GeneID	3990;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LIPC

Note: This product is for in vitro research use only